



From: Nassif, Julianne (DPH)
Sent: Friday, December 05, 2008 2:41 PM
To: Salemi, Charles (DPH); Hanchett, James (DPH)
Subject: FW: Chromatoprobe - Agilent GC & GC/MS

may be worth investigating..see Paul's note below.
Julie

From: Servizio, Paul (DPH)
Sent: Friday, December 05, 2008 2:38 PM
To: Nassif, Julianne (DPH)
Subject: FW: Chromatoprobe - Agilent GC & GC/MS

Slick trick for \$\$ 5K.. Sounds better than the old type DI probes. DLab used a probe in old days for LSD. Was Good for inlet sensitive compounds as LSD. Could be another part of the arsenal for unknown investigations. For CT.

From: OnayTech Sales [mailto:sales@onaytech.com]
Sent: Monday, December 01, 2008 12:00 PM
To: OnayTech Sales
Subject: Chromatoprobe - Agilent GC & GC/MS



The OnayTech Chromatoprobe direct sample introduction device is a revolutionary inlet for directly measuring solid samples - without any sample preparation. The Chromatoprobe easily and quickly adapts to your existing Agilent 4890, 5890, 6890 or 7890 split/splitless inlet.

The installation takes approximately 2 minutes and is as easy to install as changing the inlet liner. The Chromatoprobe can be used to quickly analyze a variety of solid samples. Some of the most popular analysis include:

- **Street drug confirmation**
 - suspicious powders
- **Natural products**
 - coffees, teas, herbals and medicinals
- **Pesticides**
 - environmental samples and food products
- **Polymers**
- **Adulterated food products**

For more information, visit us at
www.onaytech.com/cprobe.html